


THROUGH THE SKIN IN A FRACTION OF A SECOND.

With the INJEX30, the insulin is pushed through a micro-orifice into the superficial skin cells to the tissue below. This is achieved by using a precise measurement of pressure at a fraction of a second.

This is accomplished by using a spring.





Light in Weight

The INJEX30 is travel friendly. The injector which is the size of a pencil (only 96 mm), along with the reset box and accessories, fits conveniently into a small carrying case.

Long life

An enormous advantage of the INJEX30 is its long life. The injector is made of durable stainless steel and is good for approximately 5,000 injections. In addition to its small size (96 mm) and weight (75 g), it comes with one year warranty.




THE INJEX30 SYSTEM

a. Injector (The syringe without needle)

b. Transporter (dosing tool, designed for transferring liquid medication from cylinder cartridges to the INJEX system ampoules)

c. Resetbox (for charging the Injector)



INJEX30 ACCESSORIES

a. SiliTop (for dental application)

b. Ampoule (drug delivery disposable sterile medication cartridge)

c. Several Adapters (for filling the Ampoules from several devices)

0,17 mm




0,39 mm



LESS PAIN

Use of the INJEX30 needle free injection system is less painful because the skin is penetrated by an extremely fine jet of medication. This is a tremendous help to patients who are uncomfortable or fear the use of traditional needles.

The puncture point is considerably smaller [see illustration].



SAFE

Trials with this system have shown that the extremely fine jet of medication gently penetrates the subcutaneous fatty tissue and selects the path of least resistance. In effect there is no damage to blood vessels, nerve fibers or osseous tissues have been observed.



INJEX30 FOR INSULIN THERAPY

With the INJEX30, the insulin is pushed through a micro-orifice into the superficial skin cells to the tissue below. This is achieved by using a precise measurement of pressure at a fraction of a second.

0,17 mm




0,39 mm



LESS PAIN

Use of the INJEX30 needle free injection system is less painful because the skin is penetrated by an extremely fine jet of medication. This is a tremendous help to patients who are uncomfortable or fear the use of traditional needles.

The puncture point is considerably smaller [see illustration].



SAFE

Trials with this system have shown that the extremely fine jet of medication gently penetrates the subcutaneous fatty tissue and selects the path of least resistance. In effect there is no damage to blood vessels, nerve fibers or osseous tissues have been observed.



INSULIN USE DOES NOT CHANGE

The INJEX30 can be used with all licensed types of insulin.

THE REQUIRED UNITS OF INSULIN MUST NOT BE ALTERED WHEN YOU CHANGE TO INJEX.

Insulin with different periods of effect can be mixed (mixture as recommended by the manufacturer).